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catalog

| | |
|---|----|
| 1. Intended use、 Model、 Applicable strip、 Test Principle .. | 3 |
| 1.1 Intended use | 3 |
| 1.2 Model | 3 |
| 1.3 Applicable test strip | 3 |
| 1.4 Test Principle | 3 |
| 2. Structure and Diagram of Meter and Accessories | 3 |
| 2.1 Main structure composition | 3 |
| 2.2 Diagram of Meter and Accessories | 4 |
| 3. Display | 7 |
| 4. Operation | 11 |
| 4.1 Changing Battery | 11 |
| 4.2 Meter set up | 12 |
| 4.3 Blood Glucose Testing | 14 |
| 4.4 Meter Memory | 16 |
| 5. Power saving function | 18 |
| 6. Check The System | 18 |
| 6.1 Automatic Self-Test | 18 |
| 6.2 Control Test | 20 |

| | |
|---|----|
| 7. Bluetooth | 21 |
| 8. Preservation and Maintenance | 22 |
| 9. Warranty | 23 |
| 10. Matters need attention | 23 |
| 11. Contraindications | 24 |
| 12. Troubleshooting | 24 |
| 13. Performance | 29 |
| 14. Specification | 30 |
| 15. Accessory list | 32 |
| 16. Explanation of symbols | 32 |
| 17. EMC Statement | 32 |
| 17.1 Guidance and manufacturer' s declaration -electromagnetic emission | 34 |
| 17.2 Guidance and manufacturer' s declaration -electromagnetic immunity | 35 |
| 17.3 Guidance and manufacturer' s declaration -electromagnetic immunity | 36 |
| 17.4 Recommended separation distances between portable and mobile RF communications equipment and the equipment | 37 |

1.Intended use, Model, Applicable strip, Test Principle

1.1 Intended use

The product can be used with applicable blood glucose test strip to test glucose in whole blood sample. It can be used at clinical setting or self-testing by diabetes and other groups. The product is only use for monitoring glucose, not for use in diagnosis of diabetes mellitus.

1.2 Model

TRUE METRIX、 TRUE METRIX AIR.

1.3 Applicable test strip

TRUE METRIX blood glucose test strip (if need to buy test strip, contact sinocare or its authorized dealer to obtain it)

1.4 Test Principle

Blood glucose test strip reaction zone has specific chemicals, glucose in whole blood sample react with these chemicals and produce an electrical current. The meter measures the current and calculates the amount of glucose. The result is diplayed as a plasma value.

2.Structure and Diagram of Meter and Accessories

2.1 Main structure composition

Meter: Meter is composed by PCB board, buttons, LCD display and case.

Accessories: Lancing device (product which already get registration certificate or filing), batteries, docking station (optional), USB cable (optional).

Note: Docking station (optional) only apply to 真睿TRUE METRIX、真睿TRUE METRIX AIR. USB cable (optional) only apply to 真睿TRUE METRIX GO.

2.2 Diagram of Meter and Accessories

2.2.1 真睿TRUE METRIX、真睿TRUE METRIX AIR Schematic diagram

① “ ◀ ” Button

② “ • ” Button

③ “ ▶ ” Button

真睿TRUE METRIX:

真睿TRUE METRIX AIR:

① Display Screen

② Test Port

③ Strip Release Button

④ Battery Door


④ Battery Door

⑤ Meter Label

⑥ Docking Station Contacts

(It is used for data transmission with PC after connecting Docking station. Docking station is an optional accessory. If you need to buy Docking station, please contact sinocare or its authorized dealer to obtain it)

2.2.2 Lancing device


Please refer to lancing device instruction of use for more information.

Note: Blood glucose meter can only use with matched strip. Do not use with other product.

3.Display

Full screen display

1.Result is from Memory

2.Time, Date

3.Bluetooth Symbol
(only for 真睿 TRUE METRIX
AIR) and Event Tag Symbols

4.Average value symbol

5.Test Result

6.Control Symbol

7.Battery Symbol

8.Alternate Site (ALT) Symbol

9.Drop Symbol

10.Temperature Symbol

11.Test Reminder Symbol

12.Ketone Test Alert Symbol

13.Units Symbol

Waiting for apply sample

During testing

Blood glucose test result (example)

Control solution test result (example)

Event tag

Test reminder on

Test reminder off

Ketone alert on

ALT testing

Setting year


Setting date


Setting time

No average result


7-day average
test result

14-day average
test result

30-day average
test result

60-day average
test result
(Only for 真睿 TRUE METRIX AIR)

90-day average
test result
(Only for 真睿 TRUE METRIX AIR)

Blood glucose test
result in memory
(sample)

Control solution test
result in memory
(sample)

Bluetooth
function on
(Only for 真睿 TRUE METRIX AIR)

Bluetooth
function off
(Only for 真睿 TRUE METRIX AIR)

Invalid Hematocrit

Temperature error
Too cold/Too hot

Sample not detected
or used wrong
test strip

Used test strip
or apply the
sample error

Meter error

Test strip error

Test strip removed
during test.

Communication error

Low battery

Broken Display

Out of range -
High results > 33.3 mmol/L

Out of range -
Low results < 1.1 mmol/L

4. Operation

4.1 Changing Battery

 Do not recharge battery. USB cable is only for data transmission. Do not use USB cable to recharge. Battery may explode if mishandled. Do not dispose of battery in fire. Do not take apart or attempt to recharge battery. Don't throw away, dispose according to local regulations.

The meter use one 3V CR2032 lithium battery (non-recharge). Replace battery when display shows low battery error.

1. Lift tab on Battery Door.

2. Turn meter over, dispose according to local regulations.

3. Insert new battery, positive ("+") side facing up. Close Battery Door

4. Press "•" Button to turn meter on. If meter does not turn on, check that battery was installed properly. Install the battery cover and repeat step 4. Contact for assistance if problem persists.

Note: If battery is out of meter or low battery, meter may reset to original factory settings. Verify settings are correct after replacing battery by, going into Meter Set Up and checking time, date, Ketone Testing Alert, and Testing Reminders. If not, please enter the setting mode and change the setting. Results in Memory are not deleted, if battery is low or removed for any length of time.

4.2 Meter set up

Meter comes with pre-set time and date. The Event Tags, Ketone Test Alert and all Test Reminders are off. Bluetooth function is initially turned on (only for 真睿 TRUE METRIX AIR). Before using the meter for the first time or after a battery change, please check the time, date. Follow these steps to update the Settings if needed.

1. With meter off, press and hold “•” Button until the full Display is shown and a series of beeps sound. Release “•” Button. Meter goes into Set Up.

Note: If the meter turns off at any time during Set Up, go back to Step 1 and begin again.

Set Time/Date

1. The hour flash after enter the setting mode. Press “▶” or “◀” Button on top of the meter to select the hour. Press “•” Button to set.

2. The minute flash after setting hour. Press “▶” Button or “◀” to select the minutes. Press “•” Button to set.

3. Repeat step 2 to set month, date and year.

Note: Meter beeps every time a setting is confirmed.

Set Bluetooth Smart (only for 真睿 TRUE METRIX AIR)

Uploading results using the Bluetooth Smart feature allows meter to send results by connecting Bluetooth device.

After setting the year, press “▶” or “◀” Button to turn on or off the Bluetooth Smart feature. Press “•” Button to set. Meter goes to set Event Tags.

Event Tags

Event Tags are used to mark a test result that was taken during a specific event. Event tagging may assist you and your doctor or healthcare professional with managing your diabetes.

After turning on or off the Bluetooth Smart feature (for 真睿 TRUE METRIX AIR meter) or after setting the year (for 真睿 TRUE METRIX meter), press "▶" or "◀" Button to turn Event Tags on or off. Press "•" Button to set. The meter goes to set Ketone Test Alert.

-  Before meal –test was taken just before a meal;
-  After meal –test was taken 2 hours after the start of a meal;
-  Exercise – test was taken during or just after exercise;
-  Medications – test was taken after medication;
-  Sick – test was taken when sick;
-  Other – any other reason that the test is unique or different in some way (example: stress, drinking alcohol).

Ketone Test Alert

When a blood glucose result is over 13.3 mmol/L, the Ketone Test Alert is a reminder to check ketones per the treatment plan.

Press "▶" or "◀" Button to turn Alert on or off. Press "•" Button to set. Meter goes to set Test Reminder.

Note: When a Ketone Test Alert Symbol appears, it does not mean that ketones have been detected in the blood. Perform a ketone test per the treatment plan, as prescribed by the doctor or healthcare professional

Test Reminder

Up to four Test Reminders per day may be set. Meter reminder sounds at set time for 10 seconds. Test reminders are initially turned off.

1. After pressing “•” Button to set Ketone Test Alert, display shows first Reminder setting. Press “▶” or “◀” Button to turn Reminder on or off. Press “•” Button to set.

2. When “on” is chosen, press “•” Button. The hour flashes. Press “▶” or “◀” Button to set the hour. Press “•” Button to set.

3. The minutes flash. Press “▶” or “◀” Button to set the minutes. Press “•” Button to set. Meter goes to the next Test Reminder.

4. Turn Reminders on and repeat steps 1-3 to set the time for next 3 Reminders (if needed).

Note: If Test Reminders are set, the Test Reminder Symbol appears in all Displays.

Exit setting mode

Press and hold “•” Button until meter turns off. Meter also turns off after 2 minutes of non-use. Set-up choices are saved.

4.3 Blood Glucose Testing

Refer to lancing device Instructions for Use for detailed instructions.

1. Check dates on test strip vial being used. Do not use if either 4 months after first opening or after date printed next to on label, whichever comes first. And use new test strip for testing.

2. Allow vial of test strips and meter to adjust to room temperature for 10 minutes.

3. Clean hands and area to be lanced with an approved disinfectant. Dry thoroughly.

4. Remove one test strip from vial. Close vial immediately. Use test strips quickly after removal from vial.

5. With meter off, insert a test strip Contact End (blocks facing up) into Test Port. Meter turns on. Keep test strip in meter until testing is finished. To mark test as alternate site (forearm) result if needed, press "▶" Button. ALT Symbol appears in Display. Press "◀" Button to remove ALT Symbol.

6. Wait until Drop Symbol appears in Display. Do not apply samples before insert strip into meter.

7. Lance fingertip or forearm. Obtain a blood sample, allow drop to form.

8. With test strip still in meter, touch Sample Tip of test strip to top of blood drop and allow blood to be drawn into test strip.

9. Remove Sample Tip from blood drop immediately when the meter beeps.

10. After the test is finished, result is displayed. Recording test results.

11. To mark the result with an Event Tag, Event Tags must be turned on. The Event Tag icons flash. Press "▶" or "◀" Button to go to the correct Event Tag. Press "•" Button to set, event Tags stop flashing . Event tags can only be marked before the strip is removed.

12. Hold meter with test strip pointing down. Press Strip Release Button to discard test strip in the appropriate container. Meter turns off. Result is stored in Memory. The used test strip and lancet should be treated as medical waste.

Note:

-If test strip has been out of the vial too long before testing, an error message appears. Remove and discard old test strip. Use new test strip for testing.

-Do not place blood drop on top of test strip.

-Removing the test strip before result is displayed cancels the test. An error message appears. Result is not stored in Memory. Retest with a new test strip. Do not remove before result is displayed.

-Only professional can collect venous samples.

System and Laboratory Testing

When comparing results between the product and a laboratory system, meter blood tests should be performed within 30 minutes of a laboratory test.

4.4 Meter Memory

View Averages

1. With meter off press and release “ • ” Button.
2. Display scrolls through 7-, 14-, and 30-Day average values for 真睿 TRUE METRIX meter. Display scrolls through 7-, 14-, 30-, 60 and 90-Day Average values for 真睿 TRUE METRIX AIR meter. Without any operation for about 2

minutes, the blood glucose meter automatically enters the interface of historical results query. If there is no average, the glucose meter shows three dashes.

3. Meter turns off after 2 minutes if no buttons are pressed.

View Results

Meter Memory stores 500 results for 真睿 TRUE METRIX meter , When the 501th blood glucose test result is stored, the earliest blood glucose test result is automatically covered.

Meter stores 1000 results for 真睿 TRUE METRIX AIR meter. When the 1001th blood glucose test result is stored, the earliest blood glucose test result is automatically covered.

1. Press and release " • " Button after meter displays Averages value.

2. The screen shows the latest blood glucose test value. Blood glucose test results are displayed together with storage results symbol, time and date.

3. Press "▶" Button or "◀" Button to Browse test results .

Note: If test results marked as alternate site display ALT Symbol, test results display the ALT Symbol. Control Test results display the Control Symbol. If the Ketone Test Alert function on, when test results above 13.3 mmol/L, test results display Ketone Test Alert Symbol. If test results mark the event tags, test results display the event tags.

5. Power saving function

1. The waiting time for applying sample is about 2 minutes, the error of E-2 is given if over time, and the meter turns off automatically after about 10 seconds.

2. After the test result appears, if the test strip is not pulled out and there is no operation, the meter turns off automatically after about 2 minutes.

3. During the result query, there is no operation under the average result query interface, and it will automatically enter the historical result query interface after about 2 minutes. There is no operation under the historical result query interface, and it will automatically turn off after about 2 minutes.

4. When setting the blood glucose meter, if there is no operation, it will automatically turn off after about 2 minutes.

5. Automatically shut down after about 10 seconds under the condition of error reporting.

6. Check The System

Daily quality control as a means of testing the running state of the blood glucose test system. The System offers two kinds of quality control tests to let you know that the system is working properly and testing technique is good. Regular systematic inspection can ensure accurate results of blood glucose meter.

6.1 Automatic Self-Test

Each time the test strip is inserted into the blood glucose meter, the blood glucose meter starts to automatically measure itself.

Note: The Automatic Self-Test does not take the place of running a Control Test.

1. Insert a test strip into the Test Port.

2. If the blood glucose meter works normally, the blood glucose meter will display in full screen, the time will be displayed on the top of the screen, and the symbol of blood drop will blink, check for missing segments.

3. If an error message appears, do not use the blood glucose meter for testing and see Troubleshooting.

6.2 Control Test

 The quality control range affixed to the test strip bottle is not the recommended blood glucose level, and is only used for reference to the quality control solution test results. Do not take the quality control solution.

6.2.1 When to perform control test

- before using the meter for the first time;
- for practice to ensure your testing technique is good;
- when opening a new vial of test strips;
- occasionally as a vial of test strips is used;
- if results seem unusually;
- whenever a check on the glucose meter or test strips is working properly;
- if meter damage is suspected (meter was dropped, crushed, wet, etc.);

Important note: Three levels of control solution are available for Control Tests. Please use at least 2 concentrations of control solution to test whether the system is working properly.

6.2.2 How to Test Control Solution

Use ONLY 真睿 TRUE METRIX Control Solution for Control Test.

1. Check dates on control solution label and test strip vial label. Do not use control solution if control solution - 3 months after opening or date next to on label (whichever comes first), and test with a new blood glucose control solution. Do not use test strips - 4 months after opening or date next to symbol (whichever comes first), and use new test strip for testing.

2. Allow control solution, vial of test strips and meter to adjust to room temperature for 10 minutes before use. If it is the first time to open the blood glucose quality control bottle, please mark the date on the bottle label.

3. Gently swirl or invert control solution bottle to mix. DO NOT SHAKE! Remove one test strip from vial. Close test strip vial immediately. Use test strip quickly after removal from vial.

4. Insert test strip into Test Port. Meter turns on. Keep test strip in meter until testing is finished. Do not add control solution to test strip before inserting into meter.

5. Remove cap from control solution bottle. Gently squeeze a drop onto a clean tissue. Wipe off bottle tip and discard tissue. Gently squeeze a drop onto a small piece of unused aluminum foil.

6. Touch Sample Tip of test strip to top of drop of control solution. Allow drop to be drawn into test strip.

7. Remove test strip from drop when meter begins testing.

8. After testing is finished, result appears in the meter Display with the Control Symbol.

9. When the glucose meter shows results, compare result to Control Test range printed on test strip vial label for level of control solution you are using. If result is in range, system can be used for testing blood. If result does not fall within range, repeat control test using a new test strip. If result is still outside range, do not use this blood glucose meter and test strips to test. Please call the after-sales service.

10. After result is shown, Hold glucose meter with test strip pointing down, Press Strip Release Button to discard test strip, And according to medical waste disposal. Meter turns off. Recap control solution bottle tightly.

7. Bluetooth

真睿TRUE METRIX AIR glucose meter has Bluetooth communication function.

The Bluetooth function can be turned on or off by setting the blood glucose meter. When the Bluetooth function is turned on, the Bluetooth function will be turned on and the Bluetooth symbol will appear on the display interface of the glucose meter. When the Bluetooth connection is successful, the Bluetooth symbol is always on. The Bluetooth symbol continues to flash if the connection is not successful.

When the Bluetooth connection is successful, the blood glucose meter will upload the test result to Bluetooth device. Bluetooth device will display the test result when upload successfully; and the test result will be flicker at most 5 seconds before the blood glucose meter power off when upload failed

In the state of result query, if the Bluetooth connection is successful, the historical results stored by the glucose meter that have not been successfully uploaded to the Bluetooth device will be uploaded again.

When the glucose meter is turned off, the Bluetooth function will also be turned off. No device connection allowed.

8. Preservation and Maintenance

1. Store meter in carrying case to protect from liquids, dust and dirt. Do not keep meter in an area where it may be crushed (i.e. back pocket, drawer, bottom of bag, etc.).

2. Clean immediately after getting any blood on the meter or if meter is dirty or before allowing anyone else to handle it. Do not use bleach to clean the meter. Make sure meter is off and a test strip is not inserted. Wipe Meter with a clean, lint free cloth or swab dampened with 75% alcohol or neutral cleaner. Let Meter air dry thoroughly before using to test.

3. If not use meter for long time, wipe Meter with a clean, lint free cloth or swab dampened with 75% alcohol or neutral

cleaner. Let Meter air dry thoroughly and remove battery.

4. Do not clean the meter during a test.

5. Never put Meter in liquids or allow any liquids to enter the Test Port.

9. Warranty

Only Sinocare or Sinocare 's agency can inspect or supply any spare parts.

In normal use, if the product fails, the company promises 10 years warranty. Please fill in the warranty card carefully and return the parts that need to be returned to our company.

If blood glucose meter need to return to our company in case of failure, use wipe or swab with 75% alcohol to clean meter then send back to our company.

10. Matters need attention

1. Use of meter in a manner not specified in this Owner' s Booklet is not recommended and may affect ability to determine true blood glucose levels.

2. The product is an in vitro (outside body) quantitative system that is used for self-testing of human whole blood only.

3. Do not use this product during a xylose absorption test.

4. The test results of this blood glucose meter can only be used for blood glucose monitoring, and cannot be used as the basis for diabetes diagnosis. NEVER change a treatment

plan without consulting a doctor or healthcare professional. When you use this blood glucose meter to get test results inconsistent with the symptoms, you should immediately go to the hospital for examination.

5. DO NOT perform testing on the critically ill or newborn.

6. During the test, the blood glucose meter may come into contact with blood. There is therefore a risk of infection with a used glucose meter. When the glucose meter is used in a medical setting, health care workers should follow the appropriate infection control procedures for your health equipment, such as wearing gloves or other personal protection.

7. Children or special users in need of monitoring must use this blood glucose meter under the supervision of other normal adults, and the blood glucose meter must be placed out of the reach of children.

8. In order to avoid electromagnetic interference from other devices affecting the test results of this product, please keep the measurement away from strong interference sources such as short-wave high-frequency equipment.

11. Contraindications

NONE

12. Troubleshooting

1. After inserting test strip, meter does not turn on.

| Reason | Action |
|--|---|
| Test strip inserted upside down or backwards | Remove test strip. Re-insert correctly. |
| Test strip not fully inserted | Remove test strip. Re-insert test strip fully into meter |
| Test strip error | Repeat with new test strip. |
| Low battery | Replace battery. |
| Meter error | Reinstall the battery, let battery positive ("+") side face up. |
| Battery in backwards | Contact for assistance. |

2. After applying sample, test does not start.

| Reason | Action |
|--|--|
| Sample drop too small | Repeat test with new test strip and larger drop. |
| Sample applied after two minute shut-off | Repeat test with new test strip. Apply sample within 2 minutes of inserting test strip. |
| Test Strip Error | Repeat with new test strip. If you are still unable to start the test, please contact for assistance |
| Meter Error | Contact for assistance. |

3. Error messages

| Display | Reason | Action |
|--|---|---|
|  | Invalid Hematocrit | Repeat with new test strip, using capillary whole blood from the finger. If error persists, contact for assistance. |
|  | Temperature error Too cold/Too hot | Move meter and test strips to area between 5°C-40°C, wait 10 minutes for system to reach room temperature before testing. |
|  | Sample not detected or using wrong test strip | Retest with new test strip according to the instructions. |
|  | Used test strip | Repeat with new test strip. If error persists, contact for assistance. |

| Display | Reason | Action |
|--|--------------------------------|---|
|  | Meter error | Contact for assistance. |
|  | Test strip error | Retest with new test strip. If error persists, contact for assistance. |
|  | Test strip removed during test | Retest with new test strip. Make sure result is displayed before removing test strip. |
|  | Communication error | Contact for assistance. |

| Display | Reason | Action |
|--|---|---|
|  | Low battery | Change new battery. |
|  | Broken Display | Do not use meter for testing. Contact for assistance. |
|  | Out of range - High results > 33.3 mmol/L | Retest with new test strip. If result is still "Hi" (High) contact doctor or healthcare professional immediately. |
|  | Out of range - Low results < 1.1 mmol/L | Retest with new test strip. If result is still "Lo" (Low) contact doctor or healthcare professional immediately. |

13. Performance

Accuracy and repeatability criteria for this product is refer to national standard GB/T19634-2005 In vitro diagnostic test systems -General technical requirements for blood-glucose monitoring systems for self-testing and International standard EN ISO 15197:2015 In vitro diagnostic test systems —Requirements for blood-glucose monitoring systems for self-testing in managing diabetes mellitus.

Accuracy requirement: $\geq 95\%$ of test result bias shall meet table 1 requirement.

Repeatability requirement: test result shall meet table 2 requirement.

Table 1: Accuracy Requirement

| Test Range | Accept Bias |
|--|---|
| 1.1mmol/L ~ 5.55mmol/L (20mg/dL ~ 100mg/dL) | Within $\pm 0.83\text{mmol/L}$ ($\pm 15\text{mg/dL}$) |
| 5.55mmol/L ~ 33.3mmol/L (100mg/dL ~ 600mg/dL) | Within $\pm 15\%$ |

Table 2 : Repeatability Requirement

| Test Range | Precision |
|--|--------------------------------------|
| 1.1mmol/L ~ 5.55mmol/L (20mg/dL ~ 100mg/dL) | SD $< 0.34\text{mmol/L}$ (6 mg/dL) |
| 5.55mmol/L ~ 33.3mmol/L (100mg/dL ~ 600mg/dL) | CV $< 6\%$ |

14. Product performance parameters

Result Range: 1.1 - 33.3 mmol/L

Sample Size: 0.5 microlitre (0.5 μ L)

Sample: Fresh capillary whole blood, venous whole blood collected in sodium or lithium heparin blood collection tubes, or EDTA blood collection tubes, or control solution

Test Time: 4~7 seconds

Result Value: Plasma values

Assay Method: Amperometric

Power Supply: One 3V lithium battery #CR2032 (non-dis-mountable and non-rechargeable)

maximum current rating: 真睿TRUE METRIX 8 mA、真睿 TRUE METRIX AIR 10mA

Shelf life: 10 years (based on every day perform 7 tests)

The use period of this product is based on every day perform 7 tests. During the use, users should maintain and maintain the product according to the requirements of the product manual. After maintenance, products that are confirmed to maintain basic safety and effectiveness can be used normally.

Battery Life: Approximately 1000 tests or 1 year

Manufacture date: Refer to meter label

Software release: 真睿 TRUE METRIX meter: V02

真睿 TRUE METRIX AIR meter: V02

Memory Size: 500 glucose results and 1 control result for 真睿 TRUE METRIX meter

1000 glucose results and 1 control result for 真睿 TRUE METRIX AIR meter

Size: 87mm×55mm×18mm

Operating Range: Relative Humidity: 10%-80%RH
(Non-condensing), Temperature: 5°C- 40°C

Storage and transportation conditions: 10%-93%RH
(Non-condensing), -20°C-55°C

Hematocrit: 20%-70%

Altitude: Up to 3108 meters

Network security:

1. Data interface

Serial port communication is adopted between the TRUE METRIX AIR glucose meter and the data transmission module. There are 8 data bits, 1 stop bit, no check, no flow control, and the baud rate is fixed at 9600bps. Data is transmitted in hexadecimal format.

Built-in Bluetooth (4.2) is a pass-through module designed based on Bluetooth Low Energy standard. After connection with other Bluetooth devices, two-way data transmission can be realized.

Data transmission conditions: the blood glucose meter and Bluetooth device are enabled.

2. User access control

Only by obtaining specific communication protocols can users analyze data accurately, and data exchange between glucose meters and the outside world must follow the requirements of communication protocols.

3. Operating environment

See operating instructions for software environment.

15. Accessory list

Lancing device: 1 pcs, product which already get registration certificate. Operation method please refer to lancing device instruction of use.

Battery: One 3V lithium battery #CR2032 (non-rechargeable). Details please check "4.1 Changing Battery" .

Docking station (Optional)

16. Explanation of symbols

| | |
|---|--|
|  For in vitro Diagnostic Testing Only |  Attention! Read Instructions for Use |
|  Serial Number |  Keep away from direct sunlight |
|  Keep Dry |  Electronic and electrical hazardous substance symbol |
|  Fragile, handle with care |  Manufactured By |
|  Caution, Refer to the attached document |  Biological Harm |
| CMIIT ID: Wireless meter model approval code | |

17. EMC Statement

This meter meets the electromagnetic emissions requirements as per GB/T 18268.1-2010 and GB/T 18268.26-2010. EMC basic performance: "the display screen is normally displayed, the equipment is running without fault continuously, and data transmission is normal."

Note:

(1) In dry environment, especially with artificial material (artificial textile. carpet), using this meter may cause electrostatic discharge and lead to incorrect result.

(2) Portable device or mobile radio emission device may effect performance of this meter.

(3) Do not use this meter near substantial radiation source (such as non-shielded RF sources), otherwise may effect performance of this meter.

(4) User should make sure the meter can work normally during electromagnetic compatibility environment. It is recommended to evaluate the electromagnetic environment before using the meter.

17.1 Guidance and manufacturer' s declaration- electromagnetic emission

| Guidance and manufacturer' s declaration- electromagnetic emission | | |
|---|----------------|--|
| The equipment is intended for sue in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment. | | |
| Emissions test | Compliance | Electromagnetic environment-guidance |
| RF emissions GB4824 | Group 1 | The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions GB 4824 | Class B | |
| Harmonic emissions GB 17625.1 | Not applicable | / |
| Voltage fluctuations/ flicker emissions GB 17625.2 | Not applicable | |

17.2 Guidance and manufacturer' s declaration- electromagnetic immunity

| Guidance and manufacturer' s declaration- electromagnetic immunity | | | |
|---|---|---|--|
| The equipment is intended for sue in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment. | | | |
| Immunity test | Test level | Compliance level | Electromagnetic environment-guidance |
| Electrostatic discharge (ESD) GB/T 17626.2 | Air : ±2kV;±4 kV; ±8 kV Contact : ±2 kV;±4 kV | Air : ±2kV;±4 kV; ±8 kV Contact : ±2 kV;±4 kV | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30% |
| Power frequency magnetic field GB/T17626.8 | 3A/m , 50Hz | 3A/m , 50Hz | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| Voltage dips GB/T17626.11 | 0% for cycle 1 | N/A | / |
| | 40% for cycle 5 | N/A | |
| | 70% for cycle 25 | N/A | |
| Voltage interruption GB/T17626.11 | 5%, Duration: 250 cycles | N/A | |
| Pulse train GB/T17626.4 | 1kV (5/50ns , 50kHz) | N/A | / |
| Surge GB/T17626.5 | 2kV lines to earth 1kV lines to lines | N/A | |

17.3 Guidance and manufacturer's declaration- electromagnetic immunity

Guidance and manufacturer's declaration- electromagnetic immunity

The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.

| Immunity test | Test level | Compliance level | Electromagnetic environment-guidance |
|-----------------------------|---|---------------------------|--|
| Conducted RF GB/T17626.6 | 3V (Valid values) 150kHz ~ 80MHz | N/A | Portable and mobile RF communications equipment should be used to no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2 \sqrt{P}$ 80MHz ~ 800MHz $d=2.3 \sqrt{P}$ 800MHz ~ 2.0GHz Where: P —maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer. d —the recommended separation distance in metres. |
| Radiated RF GB/T17626.3 | 3V/m 80MHz ~ 2.0GHz | 3V/m 80MHz ~ 2.0GHz | Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:  |

Note 1 : At 80MHz and 800MHz, the higher frequency range applies.

Note 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strength from fixed transmitters, such as base stations for radio telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.

b. Over the frequency range 150kHz ~ 80MHz, field strengths should be less than 3 V/m

17.4 Recommended separation distances between portable and mobile RF communications equipment and the equipment

Recommended separation distances between portable and mobile RF communications equipment and the equipment

The equipment is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the equipment as recommended below, according to maximum output power of the communication equipment.

| Rated maximum output power of transmitter W | Separation distance according to frequency of transmitter/m | | |
|---|---|------------------------------------|-------------------------------------|
| | 150kHz ~ 80MHz $d=1.2 \sqrt{P}$ | 80MHz ~ 800MHz $d=1.2 \sqrt{P}$ | 800MHz ~ 2.0GHz $d=2.3 \sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 3.8 |
| 100 | 12 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts according to the transmitter manufacturer.

Note 1 : At 80MHz and 800MHz, the higher frequency range applies.

Note 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

China type 2 diabetes blood glucose control range

(《China type 2 diabetes prevention guideline》 2017)

| State | Range |
|---|---|
| Preprandial plasma glucose (before meal) | 4.4 ~ 7.0mmol/L (79 mg/dL ~ 126 mg/dL) |

Gestational diabetes mellitus (GDM) glucose control range

(Gestational diabetes mellitus prevention guideline (2014))

| State | Range |
|---|-------------------------------------|
| Preprandial plasma glucose (before meal) | ≤ 5.3 mmol/L(≤ 95 mg/dL) |

Normal glucose range

(《The National Clinical Test Regulation of Operation 》 the fourth edition)

| State | Normal range |
|--------------|--|
| Before meal | 3.9 ~ 6.1 mmol/L (70 mg/dL ~ 110 mg/dL) |

Hazardous substance and content in product

| Component | Hazardous substance | | |
|---|--------------------------------------|--------------------------------|---------------------------------------|
| | Lead(Pb)and its compound | Mercury(Hg)and its compound | Cadmium(Cd)and its compound |
| Meter case | ○ | ○ | ○ |
| Hardware | ○ | ○ | ○ |
| PCB board | ○ | ○ | ○ |
| Electronic components | × | ○ | ○ |
| / | Hexavalent chromium and its compound | polybrominated biphenyls (PBB) | polybrominated diphenyl ethers (PBDE) |
| Meter case | ○ | ○ | ○ |
| Hardware | ○ | ○ | ○ |
| PCB board | ○ | ○ | ○ |
| Electronic components | ○ | ○ | ○ |
| <p>This table is refers to SJ/T 11364 Electrical product hazardous substance limited use regulation. ○ : Means this hazardous substance in all component meet GB/T 26572 requirement. × : Means this hazardous substance in some component do not meet GB/T 26572 requirement. (Electronic components include ceramic resistor, National suggest contamination control attestation limit substance application exception do not require limit value for ceramic resistor.) * : some products content this accessories.</p> <p> When user use this product refer to user manual, the hazardous substances in this product will not reveal or cause severe pollution to environment in 10 years.</p> | | | |

Hazardous substance and content in product

| Component | Hazardous substance | | |
|------------------------------|--------------------------------------|--------------------------------|---------------------------------------|
| | Lead(Pb)and its compound | Mercury(Hg)and its compound | Cadmium(Cd)and its compound |
| *power line or power adapter | ○ | ○ | ○ |
| *battery | ○ | ○ | ○ |
| *test strips | ○ | ○ | ○ |
| * Lancing device | ○ | ○ | ○ |
| Packaging and printing | ○ | ○ | ○ |
| / | Hexavalent chromium and its compound | polybrominated biphenyls (PBB) | polybrominated diphenyl ethers (PBDE) |
| *power line or power adapter | ○ | ○ | ○ |
| *battery | × | ○ | ○ |
| *test strips | ○ | ○ | ○ |
| * Lancing device | ○ | ○ | ○ |
| Packaging and printing | ○ | ○ | ○ |

This table is refers to SJ/T 11364 Electrical product hazardous substance limited use regulation.

○ : Means this hazardous substance in all component meet GB/T 26572 requirement.

× : Means this hazardous substance in some component do not meet GB/T 26572 requirement.

(Electronic components include ceramic resistor, National suggest contamination control attestation limit substance application exception do not require limit value for ceramic resistor.)

* : some products content this accessories.

 When user use this product refer to user manual, the hazardous substances in this product will not reveal or cause severe pollution to environment in 10 years.