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Contents

1. Intended use、 Model、 Applicable strip、 Test Principle ..	3
1.1 Intended use	3
1.2 Model	3
1.3 Applicable test strip	3
1.4 Test Principle	3
2. Structure and Diagram of Meter and Accessories	4
2.1 Main structure composition	4
2.2 Diagram of Meter and Accessories	4
2.2.1 Meter	4
2.2.2 Lancing device	5
3.Display	5
4.Operation	8
4.1 Changing Battery	8
4.2 Meter set up	9
4.3 Blood Glucose Testing	11
4.4 Meter Memory	12
5.Quality Control Testing	13
5.1 Automatic Self-Test	13

5.2 Control Test	14
5.2.1 When to perform control test	14
5.2.2 How to Test Control Solution	15
6.Preservation and Maintenance	16
7.Warranty	17
8.Matters need attention	17
9.Troubleshooting	18
10.Performance	23
11.Specification	24
12.Accessory lis	25
13.Explanation of symbols	25
14.EMC Statement	25

1. Intended use, Model, Applicable strip, Test Principle

1.1. Intended use

The meter is used with test strip to quantitatively measure glucose (sugar) in fresh capillary whole blood samples drawn from the fingertip. The System may not be used for neonates. The System is intended for at-home use outside the body (self-testing) and for use by healthcare professionals in clinical settings in order to assist in the management of diabetes. This System is not for use in diagnosis or screening of diabetes mellitus.

1.2 Model

真睿TRUE METRIX GO

1.3 Applicable test strip

TRUE METRIX blood glucose test strip

1.4 Test Principle

The Test Strip is a plastic strip containing chemicals and electrodes. When inserted into meters, glucose is measured using amperometric technology employing a glucose dehydrogenase-FAD reaction. When whole blood or control solution is drawn into the Sample Tip of the test strip, glucose in the sample reacts with the chemicals and produces an electrical current. The meter measures the current and calculates the amount of glucose. The result is displayed as a plasma value.

2. Structure and Diagram of Meter and Accessories

2.1 Main structure composition



Meter: Meter is composed by PCB board, buttons, LCD display and case.

Accessories: Lancing device (product which already get registration certificate or filing), batteries, docking station (optional), USB cable (optional).

Note: Docking station (optional) only apply to 真睿TRUE METRIX、真睿TRUE METRIX AIR. USB cable (optional) only apply to 真睿TRUE METRIX GO.

2.2 Diagram of Meter and Accessories

2.2.1 Meter

① Display

② Test Port

③ Set Button


④ Battery Tray

⑤ Meter Label

⑥ Micro USB Port

⑦ Vial Lip Cover

2.2.2 Lancing device


Please refer to lancing device instruction of use for more information.






Note: Blood glucose meter can only be used with matched strip. Do not use with other product.











3. Display

Full screen display

- ① Time, Date, Control Symbol (-C-), Alternate Site Symbol (-A-), Average Symbol (7-, 14-, or 30-day)
- ② Test Result
- ③ Memory Result
- ④ Battery Symbol
- ⑤ Units of measure - Factory set to mmol/L or mg/dL, cannot be changed by user.
- ⑥ Drop Symbol

	<p>Waiting for apply sample</p>		<p>During testing</p>
	<p>Blood glucose test result (example)</p>		<p>Control solution test result (example)</p>
	<p>Setting year</p>		<p>Setting date</p>
	<p>Setting time</p>		<p>7-day average test result</p>
	<p>14-day average test result</p>		<p>30-day average test result</p>
	<p>Blood glucose test result in memory (sample)</p>		<p>30-day average test result</p>

	<p>Temperature error Too cold/ Too hot</p>		<p>Sample not detected or Sample Drop on Top of Test Strip</p>
	<p>Used test strip; Test strip outside of vial too long</p>		<p>Sample not detected or Sample Drop on Top of Test Strip</p>
	<p>Test strip Error</p>		<p>Test strip removed during test or Micro USB Cable connected while testing.</p>
	<p>Meter error</p>		<p>Meter error</p>
	<p>Low or dead battery</p>		<p>Broken Display</p>
	<p>WARNING! Out of range - High results > 33.3 mmol/L</p>		<p>WARNING! Out of range - Low results < 1.1 mmol/L</p>

4. Operation

4.1 Changing Battery

 Do not recharge battery. USB cable is only for data transmission. Do not use USB cable to recharge. Battery may explode if mishandled. Do not dispose of battery in fire. Do not take apart or attempt to recharge battery. Dispose according to local regulations.

The meter use CR2032 lithium battery (non-recharge). Replace battery when display shows low battery error.

1. Remove meter from top of test strip vial by holding the vial and twisting the meter $\frac{1}{4}$ turn counterclockwise. Lift meter from vial top.

2. Turn meter over until the meter label is facing up. Pull Battery Tray out until battery is exposed.

3. Holding the Battery Tray over your hand, press on edge of battery until battery drops out.


4. Insert new battery into Battery Tray with "+" side facing up. Slide Battery Tray back into meter.

5. Turn meter back over and press Set Button to turn meter on. If meter does not turn on, open Battery Tray and check that the battery was inserted with the "+" side facing up. Close Battery Tray and repeat Step 5. If meter still does not turn on, call the number on the cover of this Booklet for assistance.

4.2 Meter set up

To Attach/Remove Meter to Test Strip Vial

The meter may also be used for testing without attaching to the vial.

<p>To attach:</p> <ol style="list-style-type: none">1. Set test strip vial on flat surface with vial lip facing to the left.2. With Test Port facing front, place bottom of meter firmly on vial top. Meter must be seated flat on top of vial cap.3. Holding the vial, twist the meter 1/4 turn clockwise. The Test Port area on the meter should cover the vial lip if attached properly.	
<p>To remove:</p> <ol style="list-style-type: none">1. Holding the vial, twist the meter 1/4 turn counter clockwise.2. Lift off meter off the vial top.	

Set Time/Date

The meter comes with pre-set time and date. Before using the meter for the first time or after a battery change, check time and date and update as needed.

1. Start with the meter off. Press and hold the Set Button until the full Display appears and begins to blink. Release Set Button.

2. The time appears and the hour begins flashing. Change the number by pressing the Set Button until the desired number appears.

3. After the correct number appears, the number flashes 8 times before going to the minutes.

4. Repeat steps 2-4 to set up minutes, month, day and year.

If meter turns off at any time during Set Up, go back to Step #1 and begin again.

Pressing and releasing the Set Button only makes the numbers increase by one. Once the number reaches its limit, it resets to the lowest number. Pressing and holding the Set Button scrolls the numbers. Release Set Button when desired number is reached

The meter turns on when:

- a test strip is inserted into the Test Port, or
- when Set Button is pressed and released (see Meter Memory and Time/Date Set Up).

Meter turns off when:

- the test strip is released from the meter,

- the Set Button is pressed and held for 3 seconds, or
- after 2 minutes of non-use.

4.3 Blood Glucose Testing

Refer to lancing device Instructions for Use for detailed instructions.

1. Check dates on test strip vial being used. Do not use if either 4 months after first opening or after date printed next to on label, whichever comes first.

2. Allow vial of test strips and meter to adjust to room temperature for 10 minutes.

3. Wash hands in warm, soapy water. Rinse well and dry thoroughly.

4. Remove one test strip from vial. Close test strip vial immediately. Use test strip quickly after taking it out of the vial.

5. Insert test strip firmly into Test Port. Meter turns on.

6. Wait until Drop Symbol appears in the Display. Keep test strip in meter until testing is finished. Do not add blood to test strip before inserting into meter.

7. Set lancing device aside. Lance fingertip. To help blood drop form, lower hand to waist level and gently massage from palm to fingertip. Allow blood drop to form before attempting to apply to the test strip.

8. With test strip in Meter, touch Sample Tip to top of blood drop. Allow blood to be drawn into the test strip.

9. Remove test strip from drop when dashes appear across the meter Display. Meter is testing.

10. After testing is finished, result appears in the meter Display.

11. Remove test strip from meter and discard into an appropriate waste container. Meter turns off. Result is stored in the Memory.

Note:

-If test strip has been out of the vial too long before testing, an error message appears. Remove and discard old test strip. Use new test strip for testing.

-Do not place blood drop on top of test strip.

-Removing the test strip before result is displayed cancels the test. An error message appears. Result is not stored in Memory. Retest with a new test strip. Do not remove before result is displayed.

-Only professional can collect venous samples.

System and Laboratory Testing

When comparing results between the product and a laboratory system, meter blood tests should be performed within 30 minutes of a laboratory test.

4.4 Meter Memory

View Averages (7-, 14-, and 30-Day)

1. With meter off, press and release Set Button.

2. Display scrolls through the 7-, 14-, and 30-day Averages.

Meter turns off after 2 minutes if Set Button is not pressed. If there are no Average values, three dashes are displayed for 7-, 14-, and 30-Day Averages.

View Results

Meter Memory stores 500 results which are displayed from most recent to oldest. Once Memory is full, the oldest result is replaced with the newest result.

1. After meter displays the Averages, press and release Set Button again.

2. The most recent result is shown with the Memory Symbol. Blood test results are shown with the Memory Symbol, the time and date.

3. Continue to press and release the Set Button to scroll through results. Blood test results are shown with the Memory Symbol, the time and date

Note: One Control Test result is stored in Memory. The Control Test result will show the Control Symbol in the lower left corner of the Display.

5.Quality Control Testing

To assure accurate and reliable results, the System offers two kinds of quality control tests, an Automatic Self-Test and a Control Test. These tests ensure that the system is working properly and testing technique is good.

5.1 Automatic Self-Test

Note: The Automatic Self-Test does not take the place of running a Control Test.

1. Insert a test strip into the Test Port.
2. Full Display appears. Check for missing segments.

3. Drop Symbol begins to blink. Meter may be used for testing.

If an error message appears, the meter will not perform a test. See Troubleshooting or contact for assistance.

5.2 Control Test

5.2.1 When to perform control test:

- before using the meter for the first time,
- for practice to ensure your testing technique is good,
- when opening a new vial of test strips,
- occasionally as a vial of test strips is used,
- if results seem unusually high or low,
- if the test strip vial has been left opened, exposed to extreme heat, cold, or humidity,
- whenever a check on the performance of the system is needed,
- if meter damage is suspected (meter was dropped, crushed, wet, etc.).

Note: It is important to perform Control Tests with more than one level of TRUE METRIX Control Solution. Three levels of control solution are available for Control Tests. Use contact on the cover' for more information on how to obtain control solution.

 Ranges printed on test strip vial label are for Control Test results only and are not suggested levels for blood glucose. Do not drink control solution.

5.2.2 How to Test Control Solution

Use ONLY TRUE METRIX Control Solution for Control Test.

1. Check dates on control solution label and test strip vial label. Do not use control solution or test strips if either Use By Dates have passed. (control solution - 3 months after opening or date next to on label, whichever comes first, test strips - 4 months after opening or date next to symbol, whichever comes first.) Discard expired products and use new products.

2. Allow control solution, vial of test strips and meter to adjust to room temperature for 10 minutes.

3. Wash hands, dry thoroughly.

4. Gently swirl or invert control solution bottle to mix. DO NOT SHAKE!

5. Remove one test strip from vial. Close test strip vial immediately. Use test strip quickly after removal from vial.

6. Insert test strip into Test Port. Meter turns on. Keep test strip in meter until testing is finished. Do not add control solution to test strip before inserting into meter.

7. Remove cap from control solution bottle. Gently squeeze a drop onto a clean tissue. Wipe off bottle tip and discard tissue. Gently squeeze a drop onto a small piece of unused aluminum foil, clear plastic wrap, or waxed paper for testing.

8. With test strip still in meter, touch Sample Tip of test strip to top of drop of control solution. Allow drop to be drawn into test strip.

9. Remove test strip from drop when meter beeps and begins testing.

10. After testing is finished, result appears in the meter Display with the Control Symbol.

Compare result to Control Test range printed on test strip vial label for level of control solution you are using. If result is in range, system can be used for testing blood. If result does not fall within range, repeat control test using a new test strip. If Control Test result is outside range, test again. If result is still outside range, system should not be used for testing blood. Contact for assistance (see Booklet cover for contact information).

11. After result is shown, remove test strip from meter and discard. Meter turns off. Recap control solution bottle tightly.

6. Preservation and Maintenance

1. Store meter in carrying case to protect from liquids, dust and dirt. Do not keep meter in an area where it may be crushed (i.e. back pocket, drawer, bottom of bag, etc.).

2. Clean immediately after getting any blood on the meter or if meter is dirty or before allowing anyone else to handle it. Do not use bleach to clean the meter. Make sure meter is off and a test strip is not inserted. Wipe Meter with a clean, lint free cloth or swab dampened with 75% alcohol or neutral cleaner. Let Meter air dry thoroughly before using to test.

3. If not use meter for long time, wipe Meter with a clean, lint free cloth or swab dampened with 75% alcohol or neutral cleaner. Let Meter air dry thoroughly and remove battery.

4. Do not clean the meter during a test.

5. Never put Meter in liquids or allow any liquids to enter the Test Port.

7. Warranty

Sinocare Inc. warrants this meter to be free of defects in materials and workmanship for 10 years from the date of purchase. If during this 10 year period, the meter does not work properly, Sinocare will replace it with a new meter or equivalent product free of charge. Failure of the meter due to abuse or use not in accordance with the instructions for use is not covered by this Warranty. Use wipe or swab with 75% alcohol to clean meter before send back to Sinocare.

8. Matters need attention

1. Use of meter in a manner not specified in this Owner' s Booklet is not recommended and may affect ability to determine true blood glucose levels.

2. The product is an in vitro (outside body) quantitative system that is used for self-testing of human whole blood only.

3. Do not use this product during a xylose absorption test.

4. NEVER change a treatment plan without consulting a doctor or healthcare professional. Do not use for the diagnosis of or screening for diabetes mellitus or for measuring blood glucose in newborns. If test result do not match the way you feel, repeat the test. If the results still do not match the way you feel, call the doctor or healthcare professional.

5. DO NOT perform testing on the critically ill or newborn.

6. ALL parts of the product could carry blood-borne pathogens after use, even after cleaning and disinfecting. Reuse of devices labeled for single-use may result in product contamination and patient infection. If the meter is being operated by a second person who provides testing assistance, the meter and lancing device should be cleaned prior to use by the second person.

9. Troubleshooting

1. Use of meter in a manner not specified in this Owner's Booklet is not recommended and may affect ability to determine true blood glucose levels.

2. The product is an in vitro (outside body) quantitative system that is used for self-testing of human whole blood only.

3. Do not use this product during a xylose absorption test.


1. After inserting test strip, meter does not turn on.

Reason	Action
Test strip inserted upside down or backwards	Remove test strip. Re-insert correctly into the meter.
Test strip not fully inserted	Remove test strip. Re-insert test strip fully into meter
Test strip error	Remove test strip from meter. Repeat with new test strip.
Meter is dead or there is no battery	Remove test strip from meter. Replace battery. Use new test strip for testing.
Battery in the meter backwards	Battery must be placed in meter with positive ("+") side must face up.
Meter error	Contact for assistance.



2. After applying sample, test does not start/meter does not begin testing.

Reason	Action
Sample drop too small	Repeat test with new test strip and larger drop.
Sample applied after two minute shut-off	Repeat test with new test strip. Apply sample within 2 minutes of inserting test strip.
Test Strip Error	Repeat with new test strip. If testing still has not begun, contact for assistance.
Meter Error	Contact for assistance.

3. Messages

Display	Reason	Action
	Invalid Hematocrit	Repeat with new test strip, using capillary whole blood from the finger, forearm or venous whole blood collected only in a sodium heparin blood collection tube. If error persists, contact for assistance.
	Temperature error Too cold/Too hot	Move meter and test strips to area between 5°C-40°C; wait 10 minutes for system to reach room temperature before testing.
	Sample not detected or Sample Drop on Top of Test Strip	Retest with new test strip and larger sample. Make sure Sample Tip of test strip touched top of sample drop.
	Used test strip; Test strip outside of vial too long.	Repeat with new test strip. If error persists, contact for assistance.

Display	Reason	Action
	Meter error	Contact for assistance.
	Test strip Error or Very High Blood Glucose Result (higher than 33.3 mmol/L)	Retest with new test strip. If error persists, contact for assistance. If you have symptoms such as fatigue, excess urination, thirst or blurry vision, follow a doctor or healthcare professional's advice for high blood glucose.
	Test strip removed during test or Micro USB Cable connected while testing.	Unplug Micro USB cable. Retest with new test strip. Make sure result is displayed before removing test strip. If error persists, contact for assistance.
	Meter error	Contact for assistance.

Display	Reason	Action
	Meter error	Contact for assistance.
	Low or dead battery	Replace battery
	Broken Display	WARNING! Retest with new test strip. If result is still "Hi" (High) contact doctor or healthcare professional immediately.
	WARNING! Out of range - High results > 33.3 mmol/L	Repeat with new test strip. If error persists, contact for assistance.

Display	Reason	Action
	<p>WARNING! Out of range - Low results < 1.1 mmol/L</p>	<p>WARNING! Retest with new test strip. If result is still "Lo" (Low) contact doctor or healthcare professional immediately.</p>

10. Performance

Accuracy and repeatability criteria for this product refer to national standard GB/T19634-2005 In vitro diagnostic test systems -General technical requirements for blood-glucose monitoring systems for self-testing and International standard EN ISO 15197:2015 In vitro diagnostic test systems —Requirements for blood-glucose monitoring systems for self-testing in managing diabetes mellitus.

Accuracy requirement: $\geq 95\%$ of test result bias shall meet table 1 requirement.

Repeatability requirement: test result shall meet table 2 requirement.

Table 1: Accuracy Requirement

Test Range	Accept Bias
1.10mmol/L ~ 5.55mmol/L (20mg/dL ~ 100mg/dL)	Within ± 0.83 mmol/L (± 15 mg/dL)
5.55mmol/L ~ 33.30mmol/L (100mg/dL ~ 600mg/dL)	Within $\pm 15\%$

Table 2 : Repeatability Requirement

Test Range	Precision
1.10mmol/L ~ 5.55mmol/L (20mg/dL ~ 100mg/dL)	SD < 0.34mmol/L (6 mg/dL)
5.55mmol/L ~ 33.30mmol/L (100mg/dL ~ 600mg/dL)	CV < 6%

11. Specification

Result Range: 1.1 - 33.3 mmol/L

Sample Size: 0.5 microliter (0.5 μ L)

Sample: Fresh capillary whole blood or control solution

Test Time: Results in as fast as 4 seconds

Result Value: Plasma values

Assay Method: Electrochemical

Power Supply: One 3V lithium battery #CR2032 (non-re-chargeable), Maximum rated current: 10mA

Battery Life: Approximately 1000 tests or 1 year

Manufacture date: Refer to meter label.

Shelf life: 10 years (based on every day perform 7 tests)

Software release version: V01

Memory Size: 500 glucose results

Size: 40mm \times 36mm \times 22mm

Operating Range

Relative Humidity: 10%-80% (Non-condensing)

Temperature: 5°C - 40°C

Hematocrit: 20%-70%

Storage condition: 10~93%RH, -20~55°C

Altitude: Up to and including 3108 meters

12. Accessory list

Lancing device: 1 pcs, product which already get registration certificate. Operation method please refer to lancing device instruction of use.

Battery: One 3V lithium battery #CR2032 (non-rechargeable). Details please check "4.1 Changing Battery"

USB cable : Optional

13. Explanation of symbols

 For in vitro Diagnostic Testing Only	 Attention! Read Instructions for Use
 Serial Number	 Keep away from direct sunlight
 Keep Dry	 Electronic and electrical hazardous substance symbol
 Fragile, handle with care	 Manufactured By
 Caution!	 Biological Harm

14. EMC Statement

This meter meets the electromagnetic emissions requirements as per GB/T 18268.1-2010 and GB/T 18268.26-2010.

Note:

(1) In dry environment, especially with artificial material (artificial textile. carpet), using this meter may cause electrostatic discharge and lead to incorrect result.

(2) Portable device or mobile radio emission device may effect performance of this meter.

(3) Do not use this meter near substantial radiation source, otherwise may effect performance of this meter.

(4) User should make sure the meter can work normally during electromagnetic compatibility environment.

14.1 Guidance and manufacturer' s declaration-electromagnetic emission

Guidance and manufacturer' s declaration-electromagnetic emission		
The equipment is intended for sue in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions GB4824	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions GB 4824	Class B	
Harmonic emissions GB 17625.1	Not applicable	/
Voltage fluctuations/ flicker emissions GB 17625.2	Not applicable	

14.2 Guidance and manufacturer' s declaration- electromagnetic immunity

Guidance and manufacturer' s declaration- electromagnetic immunity			
The equipment is intended for sue in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	Test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) GB/T 17626.2	Air : $\pm 2\text{kV}; \pm 4 \text{ kV};$ $\pm 8 \text{ kV}$ Contact : $\pm 2 \text{ kV}; \pm 4 \text{ kV}$	Air : $\pm 2\text{kV}; \pm 4 \text{ kV};$ $\pm 8 \text{ kV}$ Contact : $\pm 2 \text{ kV}; \pm 4 \text{ kV}$	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Power frequency magnetic field GB/T17626.8	3A/m , 50Hz	3A/m , 50Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips GB/T17626.11	0% for cycle 1	N/A	/
	40% for cycle 5	N/A	
	70% for cycle 25	N/A	
Surge GB/T17626.5	$\pm 1\text{kV}$ lines to earth $\pm 1\text{Kv}$ lines to lines	N/A	/

14.3 Guidance and manufacturer' s declaration- electromagnetic immunity

Guidance and manufacturer' s declaration- electromagnetic immunity

The equipment is intended for sue in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.

Immunity test	Test level	Compliance level	Electromagnetic environment-guidance
Conducted RF GB/T17626.6	3V (Valid values) 150kHz ~ 80MHz	N/A	Portable and mobile RF communications equipment should be used to no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2 \cdot \sqrt{P}$ 80MHz ~ 800MHz $d=2.3 \cdot \sqrt{P}$ 800MHz ~ 2.0GHz Where: P —maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer. d —the recommended separation distance in metres.
Radiated RF GB/T17626.3	3V/m 80MHz ~ 2.0GHz	3V/m 80MHz ~ 2.0GHz	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

Note 1 : At 80MHz and 800MHz, the higher frequency range applies.

Note 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strength from fixed transmitters, such as base stations for radio telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.

b. Over the frequency range 150kHz ~ 80MHz, filed strengths should be less than 3 V/m

14.4 Recommended separation distances between portable and mobile RF communications equipment and the equipment

Recommended separation distances between portable and mobile RF communications equipment and the equipment			
The equipment is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the equipment as recommended below, according to maximum output power of the communication equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter/m		
	150kHz ~ 80MHz $d = 1.2 \sqrt{P}$	80MHz ~ 800MHz $d = 1.2 \sqrt{P}$	800MHz ~ 2.0GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	3.8
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts according to the transmitter manufacturer.</p> <p>Note 1 : At 80MHz and 800MHz, the higher frequency range applies.</p> <p>Note 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

China type 2 diabetes blood glucose control range

(《China type 2 diabetes prevention guideline》 2017)

State	Range
Before meal	4.4 ~ 7.0mmol/L (79 mg/dL ~ 126 mg/dL)

Gestational diabetes mellitus (GDM) glucose control range

Gestational diabetes mellitus prevention guideline (2014))

State	Range
Preprandial plasma glucose (before meal)	≤ 5.3 mmol/L(≤ 95 mg/dL)

Normal glucose range

(National clinical laboratory operating procedures Fourth Edition)

State	Normal range
Before meal	3.9 ~ 6.1 mmol/L (70 mg/dL ~ 110 mg/dL)

Hazardous substance and content in product

Component	Hazardous substance		
	Lead(Pb)and its compound	Mercury(Hg)and its compound	Cadmium(Cd)and its compound
PCB board	○	○	○
Hardware	○	○	○
Electronic components	×	○	○
Meter case	○	○	○
/	Hexavalent chromium and its compound	Polybrominated biphenyls (PBB)	polybrominated diphenyl ethers (PBDE)
PCB board	○	○	○
Hardware	○	○	○
Electronic components	○	○	○
Meter case	○	○	○

This table is refers to SJ/T 11364 Electrical product hazardous substance limited use regulation.
 ○ : Means this hazardous substance in all component meet GB/T 26572 requirement.
 × : Means this hazardous substance in some component do not meet GB/T 26572 requirement.
 (Electronic components include ceramic resistor, National suggest contamination control attestation limit substance application exception do not require limit value for ceramic resistor.)
 * : some products content this accessories.

 When user use this product refer to user manual, the hazardous substances in this product will not reveal or cause severe pollution to environment in 10 years.

Hazardous substance and content in product

Component	Hazardous substance		
	Lead(Pb)and its compound	Mercury(Hg)and its compound	Cadmium(Cd)and its compound
*power line or power adapter	○	○	○
*battery	○	○	○
*test strips	○	○	○
* Lancing device	○	○	○
Packaging and printing	○	○	○
/	Hexavalent chromium and its compound	Polybrominated biphenyls (PBB)	polybrominated diphenyl ethers (PBDE)
*power line or power adapter	○	○	○
*battery	×	○	○
*test strips	○	○	○
* Lancing device	○	○	○
Packaging and printing	○	○	○

This table is refers to SJ/T 11364 Electrical product hazardous substance limited use regulation.

○ : Means this hazardous substance in all component meet GB/T 26572 requirement.

× : Means this hazardous substance in some component do not meet GB/T 26572 requirement.

(Electronic components include ceramic resistor, National suggest contamination control attestation limit substance application exception do not require limit value for ceramic resistor.)

* : some products content this accessories.

When user use this product refer to user manual, the hazardous substances in this product will not reveal or cause severe pollution to environment in 10 years.